

MycAway™ Spray(Ready-to-use)

Product Description

MycAway™ Spray (Ready-to-use) effectively inhibits the growth and contamination of fungi (including spores), bacteria (including spores), mycoplasma, and viruses (including HIV and Hepatitis B). Its active ingredient is a quaternary ammonium derivative. It is free from mercury, formaldehyde, phenol, alcohol, and other irritants, non-toxic, and biodegradable. It does not cause any damage to general work surfaces and equipment.

This product is a ready-to-use solution, to be used directly.

Specifications

Product Name	Product Code	Specifications
MycAway™ Spray(Ready-to-use)	C230101E	500mL
Mycoplasma Spray (Ready-to-use)	C230101S	2×500mL
	C230101M	10×500mL

Instructions

1. Spray the CO2 incubator once every two weeks. There is no need to empty the incubator before spraying.
2. Spray the laminar flow hood before each use and allow it to air dry naturally.

Shipping and Storage

Store and transport at room temperature. Shelf life is 18 months.

Notes

1. Please read the entire instruction manual carefully before using this reagent.
2. Standardized operations, including the preparation of the reaction system, sample handling, and sample addition, should be followed throughout the experiment.
3. For your safety and health, wear laboratory coats and disposable gloves during operation.

