

## Recombinant Mouse Dendritic Cell Inflammatory Protein-1/CXCL3 (Mouse DCIP-1/CXCL3)

### Product Information

Product Name	Cat#	Size
Recombinant Mouse Dendritic Cell Inflammatory Protein-1/CXCL3 (Mouse DCIP-1/CXCL3)	90923ES10	10 µg
	90923ES60	100 µg
	90923ES76	500 µg

### Product Description

CXCL3, also known as DCIP1 in mouse, CINC2 in rat, and GRO $\gamma$  in humans, is belonging to the CXC chemokine family. The amino acid sequence of mouse CXCL3 is 57 % identical to human CXCL3. CXCL3 plays a role in inflammation and exert its effects on endothelial cells in an autocrine fashion. Similar to other alpha chemokines, the three GRO proteins are potent neutrophil attractants and activators. In addition, these chemokines are also active toward basophils. CXCL3 was reportedly associated with the invasion and metastasis of various malignancies. As a member of the chemokine family, CXCL3 was previously known to participate in many pathophysiological events. CXCL3 was found to be upregulated in aggressive cancer cells.

### Product Properties

<b>Synonyms</b>	chemokine (C-X-C motif) ligand 3, CINC-2, CINC-2b, CXCL3, Dcip1, Gm1960, GRO gamma, GRO3 oncogene, GRO3, GROG, GRO-gamma, Growth-regulated protein gamma, Macrophage inflammatory protein 2-beta, MGSA gamma, MIP2B, MIP-2b, MIP2-beta, SCYB3
<b>Accession</b>	Q6W5C0
<b>GeneID</b>	330122
<b>Source</b>	E.coli-derived mouse DCIP-1/CXCL3 protein, Ala28-Ser100.
<b>Molecular Weight</b>	Approximately 7.9 kDa.
<b>AA Sequence</b>	AVVASELRCQ CLNTLPRVDF ETIQSLTVTP PGPHTCTQTEV IATLKDGQEV CLNPQGPRLQ IIKKILKSG KSS
<b>Tag</b>	None
<b>Physical Appearance</b>	Sterile Filtered White lyophilized (freeze-dried) powder.
<b>Purity</b>	>97% by SDS-PAGE and HPLC analyses.
<b>Biological Activity</b>	The biological activity determined by a chemotaxis bioassay using human CXCR2 transfected human 293 cells is in a concentration range of 10-100 ng/ml. Fully biologically active when compared to standard.
<b>Endotoxin</b>	< 1.0 EU per 1µg of the protein by the LAL method.
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4. We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom.
<b>Reconstitution</b>	Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at $\leq$ -20°C. Further dilutions should be made in appropriate buffered solutions.

### Shipping and Storage

The products are shipped with ice pack and can be stored at -20°C to -80°C for 1 year.

Recommend to aliquot the protein into smaller quantities when first used and avoid repeated freeze-thaw cycles.

### **Cautions**

1. Avoid repeated freeze-thaw cycles.
2. For your safety and health, please wear lab coats and disposable gloves for operation.
3. For research use only!