

HB220701

10× Digestion buffer 3 GMP-grade

Product Information

Product name	Cat#	Size
	10667ES03	1 mL
10×Digestion buffer 3 GMP-grade	10667ES10	10 mL
	10667ES50	50 mL

Product Description

This product is an optimized 10×Digestion buffer for different restriction enzymes to ensure optimal (100%) digestion activity. The formulation of the product: 500 mM Tris-HCl, 1 M NaCl, 100 mM MgCl₂, 1 mg/mL OsrHSA, pH 7.9 at 25°C.

This product is produced in accordance with GMP process requirements and provided in a liquid form.

Contents

Contents No.	Name	Catalog No./Specification		
		10667ES03	10667ES10	10667ES50
10667	10×Digestion buffer 3 GMP-grade	1 mL	10 mL	50 mL

Shipping and Storage

The product is shipped with dry ice and can be stored at -15°C ~ -25°C for one year.

Notes:

1. For your safety and health, please wear personal protective equipment (PPE), such as laboratory coats and disposable gloves, when operating with this product.